



Asia Hepato Gene Co. Product Datasheet

華肝基因股份有限公司 website: <http://www.ah-gene.com.tw>

806 高雄市前鎮區復興四路 12 號 8F-1 Tel : 07-9553811~62 Fax : 07-9558777

Competent cell preparation buffer

Application:

1. DH5 α , JM109, HB101, BL21 Competent cell preparation
2. DH5 α (BCRC Number: 51731) as standard, some differences between strains
3. Competent cell can be kept for at least two months

Kit Components:

Cat. No	Transform efficiency	Reaction	Components	Amount	Storage
CoCell-1	$5 \times 10^7 \sim 3 \times 10^8$ cfu/ μ g	500~100	Buffer A	2ml \times 50	25 $^{\circ}$ C
			Control plasmid(3kbp)	0.1ng	-20 $^{\circ}$ C
CoCell-2	$5 \times 10^7 \sim 3 \times 10^8$ cfu/ μ g	2600~520	Buffer A powder	8g \times 2	25 $^{\circ}$ C
			Control plasmid(3kbp)	0.1ng	-20 $^{\circ}$ C

Self prepare solution

DMSO (Dimethyl sulfoxide)

LB and SOC medium

0.5mg/ml Ampicillin

1 \times TE buffer: 10mM Tris-HCl pH8.0, 1mM EDTA

Buffer preparation:

Buffer A: Add 200ml ddH₂O to dissolve 8g Buffer A powder (Cat. No. CoCell-2),

After then filtrate with 0.22 μ m sterile filter into a sterile bottle.

Control plasmid: add 1ml 1 \times TE buffer concentration 10^{-4} ng/ μ l, store at -20 $^{\circ}$ C

Competent cell preparation:

*Make sure each step keep the cell in sterile environment

1-1. Transform efficiency $\sim 3 \times 10^8$ cfu/ μ g , 2 vial

↓ Culture 1ml LB media DH5 α (BCRC Number: 51731) in 15ml sterile centrifuge tube at 37 $^{\circ}$ C , 120 rpm for 16~17 hours.

↓ Add 16 μ l overnight culture into 1.5 ml LB media in 15ml sterile centrifuge tube

↓ Culture at 37 $^{\circ}$ C , 120 rpm for 2.5 hours (more or less time will affect transform efficiency)

↓ Pour 1.5ml DH5 α into 1.5ml sterile centrifuge tube

↓ Centrifuge 10,000 rpm for 1 min then remove LB

↓ Resuspend the pellet with 0.9ml Buffer A

↓ Stand for 20 minutes

↓ Centrifuge 10,000 rpm for 1 min then remove Buffer A

↓ Resuspend the pellet with 200 μ l Buffer A + 16 μ l DMSO

↓ Stand for 10 minutes

↓ Aliquot 2 vials (108 μ l/vial)

↓ Store at -70 $^{\circ}$ C ~ -80 $^{\circ}$ C ready for use

For research use only, not for diagnostic or therapeutic use.

1-2. Transform efficiency $\sim 3 \times 10^8$ cfu/ μ g , 65 vial

- ↓ Culture 1ml LB media DH5 α (BCRC Number: 51731) in 15ml sterile centrifuge tube at 37°C, 120 rpm for 16~17 hours.
- ↓ Add 0.5ml overnight culture into 50 ml LB media in 125ml sterile flask
- ↓ Culture at 37°C, 120 rpm for 2.5 hours (more or less time will affect transform efficiency)
- ↓ Pour 50ml DH5 α into 50ml sterile centrifuge tube
- ↓ Centrifuge 5,000 rpm for 5 min then remove LB
- ↓ Resuspend the pellet with 28ml Buffer A
- ↓ Stand for 20 minutes
- ↓ Centrifuge 5,000 rpm for 5 min then remove Buffer A
- ↓ Resuspend the pellet with 6.5ml Buffer A + 520 μ l DMSO
- ↓ Stand for 10 minutes
- ↓ Aliquot 65 vials (108 μ l/vial)
- ↓ Store at -70°C~-80°C ready for use

2. Transform efficiency $> 5 \times 10^7$ cfu/ μ g , 1 \times 5 vial

5-fold dilute $\sim 3 \times 10^8$ cfu/ μ g competent cell

- ↓ Mix 400 μ l Buffer A + 32 μ l DMSO into 1 vial (108 μ l) of $\sim 3 \times 10^8$ cfu/ μ g competent cell
- ↓ Aliquot 5 vials, each vial 108 μ l
- ↓ Store at -70~-80°C ready for use

Quality control:

	Transform efficiency	Contamination Test
1	Thaw 1 vial $\sim 3 \times 10^8$ cfu/ μ g competent cell at room temperature before use	Thaw 1 vial $\sim 3 \times 10^8$ cfu/ μ g competent cell at room temperature before use
2	Mix 1 μ l control plasmid (10^{-4} ng/ μ l) into competent cell	Mix 0 μ l control plasmid
3	4°C for 5 minutes	4°C for 5 minutes
4	42°C 45 seconds	42°C 45 seconds
5	Add 60 μ l SOC Medium	Add 60 μ l SOC Medium
6	Spread 168 μ l cells onto 20-50 μ g/ml Ampicillin plate. Excessive antibiotic significantly decreased the transformation efficiency	Spread 168 μ l cells onto 20-50 μ g/ml Ampicillin plate.
7	37°C, 12~16 hours, see 40~130 colony	37°C, 12~16 hours, should see 0 colony. But several colony contaminate won't affect target colony selection.

Medium Ingredient

LB		SOC			
Tryptone	10 g/L	Tryptone	20 g/L	MgCl ₂	0.95 g/L
Yeast extract	5 g/L	Yeast Extract	5 g/L	MgSO ₄	1.2 g/L
NaCl	10 g/L	NaCl	0.5 g/L	Glucose	3.6 g/L
		KCl	0.186 g/L		